# **cryokit**

CRYOGENIC PROTECTION

NOTA INFORMATIVA E ISTRUZIONI PER L'USO  
NOTE D'INFORMATION ET INSTRUCTIONS D'UTILISATIONS  
NOTA INFORMATIVA E INSTRUCCIONES DE USO  
INFORMATION NOTE AND INSTRUCTIONS FOR USE  
INFORMATIONSBLATT UND GEBRAUCHSANWEISUNG

CRYOPLUS-2.0  
CRYOPLUS-2.1  
CRYOKIT 300-400-550  
CRYOLITE-WP-HP-HPS

RELEASE 1/2022

INFO@CRYOKIT.NET

CRYOKIT.NET

**CRYOPLUS-2.0**  
**CRYOPLUS-2.1**

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**CRYOPLUS-2.0**  
**CRYOPLUS-2.1**

**CRYOKIT300**  
**CRYOKIT400**  
**CRYOKIT550**

These gloves are made with a special composite fabric bonded to a hydrophilic membrane. The interior features another full, breathable membrane, which does not allow liquids to penetrate inside. Protection from contact cold is guaranteed by an internal lining in soft polyester fleece.

**CE CERTIFICATION**

CE certification issued by: Ricotest Srl – Notified Lab No. 0498  
Subject to specific liquid nitrogen resistance test (-195.82°Cc), by immersion in contact for 1 minute and subsequent RCT test.

N.B.: the performance level concerns the palm side, including fingers; nevertheless, the uniformity of materials and processing makes the back of the glove also protective.

The following marking can be found on the label sewn into the inside of the gloves:

These gloves are made with special composite fabrics bonded to hydrophilic membranes. The interior features another special, breathable membrane, which does not allow liquids to penetrate inside. Protection from contact cold is guaranteed by an internal lining in soft polyester fleece.

**CE CERTIFICATION**

CE certification issued by: Ricotest Srl – Notified Lab No. 0498  
Subject to specific liquid nitrogen resistance test (-195.82°Cc), by immersion in contact for 1 minute and subsequent RCT test.

N.B.: the performance level concerns the palm side, including fingers; nevertheless, the uniformity of materials and processing makes the back of the glove also protective.

The following marking can be found on the label sewn into the inside of the gloves:

	CE marking, testifying to the conformity of the gloves to the essential health and safety requirements of Regulation (EU) 2016/425
0302	Identification no. of the Notified Body performing the inspection of the PPE manufactured pursuant (MODEL C2)
CRYOKIT	Manufacturer's registered trade mark
CRYOPLUS-2.0/2.1	Glove model
10 (example)	Size
	Pictogram referring to the information note.
abcde	EN388+pictogram relating to the mechanical risks explained below: <b>2(a) 3(b) 4(c) 2(d) x(e)</b> This pictogram indicates that the gloves were designed for mechanical risks and the numbers at the foot indicate the levels obtained during laboratory tests (see explanation on the last page).
abc	EN 511 + related pictogram for protection against cold explained below: <b>(a) 2</b> = level obtained against convective cold (0.193 m <sup>2</sup> K/W) <b>(b) 2</b> = level obtained against contact cold (0.1563 m <sup>2</sup> K/W) <b>(c) 1</b> = level obtained against water penetration (no penetration)

	CE marking, testifying to the conformity of the gloves to the essential health and safety requirements of Regulation (EU) 2016/425
0498	Identification no. of the Notified Body performing the inspection of the PPE manufactured pursuant (MODEL C2)
CRYOKIT	Manufacturer's registered trade mark
CRYOKIT400-550	Glove model
10 (example)	Size
	Pictogram referring to the information note.
abcde	EN388+pictogram relating to the mechanical risks explained below: <b>2(a) 2(b) 3(c) 1(d) x(e)</b> This pictogram indicates that the gloves were designed for mechanical risks and the numbers at the foot indicate the levels obtained during laboratory tests (see explanation on the last page).
abc	EN 511 + related pictogram for protection against cold explained below: <b>(a) 2</b> = level obtained against convective cold (0.157 m <sup>2</sup> K/W) <b>(b) 2</b> = level obtained against contact cold (0.0562 m <sup>2</sup> K/W) <b>(c) 1</b> = level obtained against water penetration (no penetration)

**CRYOLITE-WP**  
**CRYOLITE-HP**

These gloves are made of soft water-repellent cowhide grain leather. The interior features another special, breathable membrane, which does not allow liquids to penetrate inside. Protection from contact cold is guaranteed by an internal lining in soft polyester fleece. **CE CERTIFICATION**  
CE certification issued by: Ricotest Srl – Notified Lab No. 0498  
Subject to specific liquid nitrogen resistance test (-195.82°Cc), by immersion in contact for 1 minute and subsequent RCT test.

N.B.: the performance level concerns the palm side, including fingers; nevertheless, the uniformity of materials and processing makes the back of the glove also protective.

The following marking can be found on the label sewn into the inside of the gloves:

	CE marking, testifying to the conformity of the gloves to the essential health and safety requirements of Regulation (EU) 2016/425
<b>0498</b>	Identification no. of the Notified Body performing the inspection of the PPE manufactured pursuant (MODEL C2)
<b>CRYOKIT</b>	Manufacturer's registered trade mark
<b>CRYOLITE</b>	Glove model
<b>10</b> (example)	Size
	Pictogram referring to the information note.
 abcde	EN388+pictogram relating to the mechanical risks explained below: <b>2(a) 1(b) 2(c) 2(d) x(e)</b> This pictogram indicates that the gloves were designed for mechanical risks and the numbers at the foot indicate the levels obtained during laboratory tests (see explanation on the last page).
 abc	EN 511 + related pictogram for protection against cold explained below: <b>(a) 1</b> = level obtained against convective cold (0.224 m <sup>2</sup> K/W) <b>(b) 3</b> = level obtained against contact cold (0.1476 m <sup>2</sup> K/W) <b>(c) 1</b> = level obtained against water penetration (no penetration)
	EN 407 + related pictogram for protection against heat explained below: X2XXXX 2 = Level obtained against contact heat: 16 seconds at 250°C.

**CRYOLITE**  
**HPS**

These gloves are made of soft, yet tough, water-repellent cowhide grain leather. The interior features another special, breathable membrane, which does not allow liquids to penetrate inside. Protection from contact cold is guaranteed by an internal lining in soft polyester fleece. **CE CERTIFICATION**  
CE certification No. 617162701/OI issued by: Ricotest Srl – Notified Lab No. 0498  
Subject to specific liquid nitrogen resistance test (-195.82°Cc), by immersion in contact for 1 minute and subsequent RCT test.

N.B.: the performance level concerns the palm side, including fingers; nevertheless, the uniformity of materials and processing makes the back of the glove also protective.

The following marking can be found on the label sewn into the inside of the gloves:

	CE marking, testifying to the conformity of the gloves to the essential health and safety requirements of Regulation (EU) 2016/425
<b>0498</b>	Identification no. of the Notified Body performing the inspection of the PPE manufactured pursuant (MODEL C2)
<b>CRYOKIT</b>	Manufacturer's registered trade mark
<b>CRYOLITE HPS</b>	Glove model
<b>10</b> (example)	Size
	Pictogram referring to the information note.
 abcde	EN388+pictogram relating to the mechanical risks explained below: <b>3(a)32(b) 3(c) 2(d) x(e)</b> This pictogram indicates that the gloves were designed for mechanical risks and the numbers at the foot indicate the levels obtained during laboratory tests (see explanation on the last page).
 abc	EN 511 + related pictogram for protection against cold explained below: <b>(a) 2</b> = level obtained against convective cold (0.216 m <sup>2</sup> K/W) <b>(b) 2</b> = level obtained against contact cold (0.1517 m <sup>2</sup> K/W) <b>(c) 1</b> = level obtained against water penetration (no penetration)
	EN 407 + related pictogram for protection against heat explained below: X2XXXX 2 = Level obtained against contact heat: 16 seconds at 250°C.

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Kora Srl – Cryokit Cryogenic