



## Product Specification Sheet

# VWR COVERALL , SPP50G, ANTISTATIC, WHITE, WITH HOOD, ZIP DOWN, ELASITIC CUFFS, SERGED SEAM, ANKELS AND WAIST

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

*VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

### Intended use

|  |
|--|
| Food & pharmaceutical industry, lab & research |
|--|

### General Sepcification

| Item               | Standard                     |
|--------------------|------------------------------|
| Material           | 50gsm SPP Antistatic         |
| Color              | White                        |
| Elastic            | Latex free, 4mm width, white |
| Zip                | Nylon                        |
| Surface Resistance | <10 <sup>10</sup> ohm        |




# Product Specification Sheet

## Material Specifications

50gsm SPP

| Item                      | Test method   | average value standard |                             |
|---------------------------|---------------|------------------------|-----------------------------|
|                           |               | average value          | tolerance for average value |
| Weight(g/m <sup>2</sup> ) | GB/T4669-1995 | 50.00                  | +/-4                        |
| Break Tensile properties  | MD            | GB/T3923.1-1997        | >85N/5cm                    |
|                           | CD            | GB/T3923.1-1997        | >50N/5cm                    |
| Elongation rate(%)        | MD            | GB/T3923.1-1997        | >30%                        |
|                           | CD            | GB/T3923.1-1997        | >30%                        |

## Dimension Specification





## Product Specification Sheet

| Item              | Unit | Specification |          |          |          |          |          | Tolerance |
|-------------------|------|---------------|----------|----------|----------|----------|----------|-----------|
|                   |      | 113-6229      | 113-6230 | 113-6231 | 113-6232 | 113-6233 | 113-6234 |           |
|                   |      | M             | L        | XL       | XXL      | XXXL     | XXXXL    |           |
| A) Chest          | cm   | 61            | 63       | 68       | 72       | 75       | 79       | ±2.0      |
| B) Sleeve Outseam | cm   | 87            | 91       | 93       | 94       | 96       | 100      | ±2.0      |
| C) Body Length    | cm   | 87            | 88       | 95       | 97       | 102      | 105      | ±2.0      |
| D) Inseam         | cm   | 77            | 80       | 82       | 80       | 85       | 87       | ±2.0      |
| E) Zipper         | cm   | 70            | 70       | 80       | 80       | 85       | 85       | ±2.0      |
| F) Wrist Relaxed  | cm   | 8             | 8        | 8        | 9        | 9        | 9        | ±1.0      |
| G) Ankle Relaxed  | cm   | 12            | 12       | 12       | 13       | 13       | 13       | ±1.0      |
| H) Full Length    | cm   | 160           | 163      | 170      | 172      | 180      | 187      | ±2.0      |

### Packaging Specification

|                         |                         |
|-------------------------|-------------------------|
| Package/Quantity        | 1pcs/PE bag, 50pcs/Case |
| Best storage Conditions | 10° - 30°C              |
|                         | 40% - 75% humidity      |
|                         | dry                     |
| Shelf Life              | 2 years                 |
| Country of Origin       | China                   |

**Disclaimer:** VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **January 2014**