



# Specification

---

1.06560.1000 Sodium perborate tetrahydrate pure

---

	Spec. Values	
Assay (manganometric)	≥ 96.0	%
Identity	passes test	
Chloride (Cl)	≤ 0.1	%
Sulfate (SO <sub>4</sub> )	≤ 1.2	%
Heavy metals (as Pb)	≤ 0.003	%
Fe (Iron)	≤ 0.002	%

Claudia Wiegand

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.